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## Effect of N-terminal amino acid sequence on the protein synthesis using a reconstituted cell-free protein synthesis system (PUREflex®)

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
### 1. Introduction

PUREflex® is based on the PURE system.

The PURE system is a reconstituted cell-free protein synthesis system, which consists of only purified factors necessary for transcription, translation and energy regeneration.

Ref; Shimizu Y. et al. (2001) *Nat. Biotechnol.*, vol. 19, p. 751

- One of the challenges;

Differences in the synthesis efficiency depending on the proteins

- Background;

- The 3rd to 5th amino acid residues in the N-terminal region are important for efficient translation.
- Large amino acids, which are encoded by AT-rich codon, facilitate early elongation reaction.

Ref; Verma M. et al. (2019) *Nat. Commun.*, Jia B. et al. (2021) *Sci. Rep.*

- Aim of this study;



Evaluation of the effect of amino acids in the N-terminal region on the synthesis efficiency at 23, 30 and 37 °C.

### 2. Protein synthesis at 23, 30, or 37 °C

#### < Information of proteins used in this study >

Protein	M.W.	N-terminal AA Seq
DHFR <i>E. coli</i> dihydrofolate reductase	18.0 kDa	MISLIAALAV...
sfGFP Superfolder GFP	26.8 kDa	MSKGEEELFTG...
Luc Firefly luciferase	60.7 kDa	MEDAKNIKKG...
ALP <i>E. coli</i> alkaline phosphatase	47.3 kDa	MRTPEMPVLE...

PUREflex®2.1 (+ GSH)  
↓ + DNA  
↓ incubation at 23, 30 or 37 °C for 24 h  
↓ SDS-PAGE and quantification

The decrease in the amount of synthesized protein at 23 and 30 °C relative to 37 °C varied depending on the protein.

### 3. Influences of the second and fourth amino acid on sfGFP synthesis


PUREflex®2.1 (+ GSH)  
↓ + sfGFP DNA  
↓ incubation at 23, 30 or 37 °C for 24 h  
↓ Measurement of fluorescence or SDS-PAGE

#### < N-terminal sequence of sfGFP >


M S K G E E L F T G A A G A A . . . .  
atg tctaaggtaagaa....  
other amino acids

#### 3-1. Synthesis of sfGFP with a different amino acid at the fourth position

##### < Fluorescence of synthesized sfGFP with a different fourth amino acid >


##### < SDS-PAGE of synthesized sfGFP >


The amount of synthesized sfGFP, especially at 23 °C, differed depending on the fourth amino acid.  
When the fourth amino acid was a large amino acid such as Tyr or His, the amount of sfGFP was increased compared to Gly.

#### 3-2. Influences of synonymous codons at the fourth position

- The influence of amino acids on the translation efficiency at 23 °C was larger than that of synonymous codons.
- In the case of Ser, there was a difference in synthesis at 23 °C between TCN and AGN codons.

#### 3-3. Synthesis of sfGFP with a different amino acid at the second position



When the fourth amino acid was Gly, the amount of sfGFP synthesized at 23 °C differed depending on the second amino acid.

When the fourth amino acid was Tyr, there was little effect of the second amino acid except for Gly.

### 4. Effect of the N-terminal additional sequence

PUREflex®2.1 (+ GSH)  
↓ + DNA  
↓ incubation at 23, 30 or 37 °C for 24 h  
↓ SDS-PAGE and quantification

Luc MEDAKNIKKG...  
ALP MRTPEMPVLE...  
DHFR MISLIAALAV...  
SKG or SKY

Addition of Ser-Lys-Tyr at N-terminal region significantly increased the amount of synthesized ALP.

### Conclusion

- The efficiency of synthesis especially at low temperature varied from protein to protein.
- The fourth amino acid of sfGFP significantly affected the synthesis efficiency at 23 °C.
- Addition of Ser-Lys-Tyr to N-terminal region significantly increased the amount of synthesized ALP at low temperatures.

- Why does the fourth amino acid affect translation at low temperature?
- Would the same result be obtained when expressed in *E. coli*?