

In vitro affinity maturation based on Ribosome Display System with PURE*frex*


Kanehisa Kojoh; Kumiko Tsuihiji; Shizue Katoh; Takashi Kanamori

GeneFrontier Corporation

Abstract


The affinity maturation is one of the key technologies to distinguish the binder from the others against the same target. For easy and feasible affinity maturation, ribosome display is well known to give great advantages. Especially, ribosome display on the PURE system was reported as a powerful method not only for *in vitro* selection of specific binders such as antibodies or designed proteins like scaffolds, but also for *in vitro* affinity maturation. The PURE system is one of the cell-free protein synthesis system reconstituted from the purified components necessary for the translation in *E. coli*. Recently, we developed the upgraded PURE system (PURE*frex*), in which the contaminants such as RNases, lipopolysaccharides and other proteins unnecessary for the translation are extremely decreased. We also developed the ribosome display system utilizing PURE*frex* (PURE*frexRD*). The selection efficiency of PURE*frexRD* is raised approximately 100-fold or 1000-fold higher than ribosome display using the original PURE system or *E.coli* S-30 extract system. In this conference last year, we reported *in vitro* selection of the binders to Erk2 protein from our own designed protein scaffold library based on FHA domain of RNF8 (human E3 ubiquitin ligase) by PURE*frexRD*. This year, we report the application of PURE*frexRD* to *in vitro* affinity maturation by off-rate selection. We prepared two randomized RNF8 libraries based on the specific binder to Erk2 by PCR-based DNA shuffling or grafting one of the loops of RNF8. After three rounds of off-rate selection over several days, we successfully improved the affinity approximately 35-fold more than the original binder in very simple steps. On the other hand, the affinity maturation of Fab will give great advantages in terms of seamless development for generating better antibody therapeutic. Therefore, we have applied the PURE*frexRD* to the affinity maturation of Fab fragments and have developed the protocols to successfully form Fab fragments on ribosome.

①PURE*frex*:Highly purified PURE system



PURE system is a reconstituted *in vitro* transcription and translation system which consists of purified 36 proteinaceous factors and *E. coli* ribosome necessary for transcription, translation and energy recycling. It also contains amino acids, NTPs and *E. coli* tRNA, so the target protein can be synthesized just by adding the template DNA to the reaction mixture. PURE*frex* is developed as “highly purified PURE system” which has been well optimized to ribosome display; Contamination of such as Lipopolysaccharide (LPS) in the system is very low.


②mRNA-ribosome-protein complex of PURE*frexRD*



Advantages

- 1)mRNA-ribosome-protein complex is highly stable in the PURE*frex* because activities of nucleases and other inhibitory factors are very low in this system.
- 2)The mRNA encoding a scaffold protein with translation arrest sequence of SecM remains stably attached to a ribosome in the PURE*frex*.
- 3)PURE*frexRD* have high selection efficiency than using normal PURE system RD.


③*In vitro* affinity maturation with PURE*frexRD*



In off-rate selection, mRNA-ribosome-protein complexes bound to antigen-immobilized streptavidin beads are washed with a large excess of free antigen to prevent the rebinding of mRNA-ribosome-protein complexes to the beads. The higher-affinity binders can be retained on the beads for a longer washing time.

④Affinity maturation of RNF8 Erk2-binder

- Feature of RNF8 Erk2-binder
- MW : 18kDa
 - structure : β sandwich
 - Loop: 4 loop
 - Kd : 2.9×10^{-8}
 - antigen: Erk2 protein




Native structure of RNF8-peptide complex (side-view & top-view)

Last year, we reported *in vitro* selection of the binders to Erk2 protein from the designed RNF8 scaffold library. The library was constructed by randomizing the 4 loops (red lines) to NNS. The Kd value of an obtained binder is 2.9×10^{-8} . This binder was used for affinity maturation.


RNF8 is a E3 ubiquitin ligase which plays a critical role in the early DNA-damage response via the ubiquitination of histone H2A and H2AX. The global fold of the FHA domain of RNF8 is an 11-stranded β -sandwich structure with the phosphopeptide-binding surface comprised of 4 loops.

⑤Construction of mutated RNF8 Erk2-binder Library

Library 1:PCR-based DNA shuffling library (ref:Zhao et al.1998. Nat.Biotechnol., 16, 258)





Library 2:Loop4 grafting library




mRNA Libraries (1.2×10^{12} molecules)
↓
In vitro translation in 20 μ L of PURE*frex*
↓ 4°C, 30 min
Mixing with Erk2 protein (0.1 nM)-immobilized streptavidin beads
in TBST 200 μ L
↓ 4°C, 120 min
Adding of free-Erk2 protein (final 1 μ M)
↓ 4°C, incubation for long times
(1st round:1 day, 2nd round:4 days, 3rd round:7 days)
Washing five times with 500 μ L of TBST
↓
Eluting mRNA from mRNA-ribosome-protein complexes with 100 μ L of 50 mM EDTA
↓
RT-PCR→Next round selection
↓
Subcloning
↓
Selecting unique clones by sequencing
↓
Screening higher affinity clones by competitive ELISA
↓
Purifying the higher affinity clones and determining IC₅₀

⑦Selection of higher affinity RNF8 Erk2-binder clones by competitive ELISA




Competitive ELISA for evaluating the affinities of unique clones. Cell extract of *E.coli* expressing each clone was pre-incubated in absence or presence of Erk2 at the indicated concentration as competitor and then allowed to bind to immobilized RNF8 protein. After washing, remaining binders were detected by HRP-conjugated anti-FLAG antibody. The relative binding (%Binding) are indicated by the ratio between the ELISA signal in the presence of the competitor and that in the absence of the competitor. The clones showed the higher competitive inhibition (red circle) were overexpressed in *E.coli* and purified with Ni column for determination of IC₅₀.

⑧Comparison of IC₅₀ of selected higher affinity clones



IC₅₀ of the each purified clone was determined, 4 clones from library 1, and 7 clones from library 2 showed more than 10-fold higher affinity than the original RNF8 Erk2-binder. The most improved clone was Lib2-1 from Library 2 (35-fold). As a whole, the higher-affinity mutants were easily obtained from Library 2 than library 1.

⑨Affinity maturation of Fab fragment with PURE*frexRD*




In almost the cases using cell-free display methods such as ribosome display, affinity maturation of antibody are only applicable to single-chain Fvs (scFvs), but cannot applied to Fab fragment. On the other hand, it was reported that L chain shuffling method is highly effective on affinity maturation of Fab fragments with *in vitro* display method such as phage display and yeast display. Therefore, we have been developed the optimum representation of Fab fragments on ribosome. In PURE*frexRD* of Fab fragments, the randomized light chain (L chain) was presented on ribosome and the heavy chain (H chain) was translated as a free-protein.

⑩Optimization of PURE*frex* for Fab synthesis (PURE*frexSS*)


Synthesis of Fab fragment using PURE*frex* with oxidized glutathione (GSSG) and DsbC protein that catalyzed the formation of correct disulfide bounds.

- PURE*frex*
- ↓ +/- GSSG,DsbC
 - ↓ +/- H chain DNA (+FLAG)
 - ↓ +/- L chain DNA
 - ↓ incubation at 37°C for 4h
 - ↓ →Ag-coated well
 - ↓ anti-FLAG/HRP
 - ↓ +TMB



Fab fragment was synthesized using PURE*frex* supplemented with 3 mM GSSG and 23 μ g/mM DsbC protein. When both of GSSG and DsbC were added in the system, functional Fab fragment with good binding activity was synthesized.

⑪Optimization of synthesis step of H chain and L chain fragment




Condition 1 :H and L chain mRNA were translated in a tube at the same time

Condition 2 :H and L chain mRNA were separately translated in each tube, after 30 min, both tubes were mixed.

Condition 3:L chain mRNA was translated in a tube, after 30 min, H chain mRNA was added to L chain expressed tube.

- mRNA
- ↓ PURE*frexSS*
 - ↓ translation (condition 1,2,3), 37°C
 - ↓ pull down (antigen ; TNF- α)
 - ↓ recovery mRNA
 - ↓ quantitative PCR



The recovered mRNA from condition 2 and 3 was increased at 5-fold and 25-fold than condition 1, respectively. These results indicate that H and L chain mRNA should be expressed in separated tube (condition 2) or in the same tube sequentially (condition 3). In translation of condition 1, because L chain mRNAs formed the L chain mRNA-ribosome-L chain protein complex, active ribosome was decreased, and it result in reduction of H chain protein. Also when H chain mRNA was independently translated, the synthesized H chain protein was easily caused to aggregation. Therefore, it seems that the recovered mRNA from condition 2 was lower than that from condition 3.

Summary

1, We successfully demonstrated the affinity maturation with PURE*frexRD*. By means of the off-rate selection from DNA shuffling and loop grafting libraries, the affinity of mutated RNF8-binder was improved ~ 35-fold without affecting the specificity.

2, We developed the optimum display method of Fab fragments on ribosome based on PURE*frexRD*. In the optimum condition (Condition 3), functional Fab fragment was presented by 25-fold more better than in the conventional condition (Condition 1).

In conclusion, PURE*frexRD* is expected to be applied to rapid and convenient affinity maturation for developing high-affinity binders of such as antibodies or scaffold proteins, especially for therapeutic and diagnostic applications.

Contact

URL:www.genefrontier.com
E-mail:kojoh@genefrontier.com
contactus@genefrontier.com