Performance of PUREfrex® 2.0 with increased productivity and its application

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PUREfrex® 1.0 • PUREfrex® 2.0

a regular kit for the synthesis of proteins without disulfide bonds

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
a supplement for the synthesis of proteins containing disulfide bonds

DnaK Mix • GroE Mix

a supplement for the synthesis of aggregate-prone proteins

1. PUREfrex®; based on the PURE system technology

The PURE system is a reconstituted cell-free protein synthesis system, which consists of only purified factors necessary for transcription, translation and energy regeneration.

Advantage

- Low level of contamination
- Easy adjustment of the reagent composition
- PCR products usable as a template DNA

Application

- High throughput preparation of proteins (including Fab, scFv, protein toxin etc.)
- Protein science research
- *in vitro* display
- Ribosome display
- mRNA display

(Ref; Shimizu Y. et al. (2001) Nat. Biotechnol., vol. 19, p. 751)

Summary

E. coli based cell-free protein synthesis systems



	Extract system	Reconstituted system		
	S30 system	PURE system (original)	PUREfrex® 1.0	PUREfrex® 2.0
Typical Yield (µg/mL)	100-1,000	10-200	10-200	20-1,000
Contamination RNase LPS	very High very High	Low High	very Low very Low	very Low very Low
Template DNA Plasmid DNA PCR product	OK NG	OK OK	OK OK	OK OK
Customization of Reagent	Difficult	Easy	Easy	Easy
Purification of His-tagged product	OK	NG	OK	OK

2. Performance of PUREfrex 2.0


PUREfrex 2.0 has 2-10 times higher productivity than PUREfrex 1.0.

Aggregate-prone proteins can be synthesized maintaining their activity.

PUREfrex® 1.0 or PUREfrex® 2.0
↓ + template DNA (PCR product)
↓ incubation at 37°C for 4 h
↓ SDS-PAGE (1 µL of reaction/lane)
↓ CBB staining

Proteins containing disulfide bonds can be synthesized in active form.

3. Application of PUREfrex 2.0; can be applicable for synthesis/screening for Fab and scFv.

The nucleotide sequence at 5' region of HC (Heavy Chain) DNA effect on the amount of product.

PUREfrex® 2.0 + DS supplement

↓ + HC DNA (PCR product) from Herceptin
↓ incubation at 37°C for 4 h
↓ SDS-PAGE (0.5 µL of reaction/lane) and Oriole staining

a) nucleotide sequence of HC DNA



	M	Q	V	L	V	S	G	C	(%)
E. coli freq. codon	atg	cag	gtg	cag	ctg	gtg	gaa	agc	ggc
	gaa								59%
Q-1	atg	cag	gtg	cag	ctg	gtg	gaa	tgg	ggc
Q-2	atg	cag	gtg	cag	ctg	gtg	gaa	tgg	ggc
Q-3	atg	cag	gtg	cag	ctg	gtg	gaa	tgg	ggc
Q-4	atg	cag	gtg	cag	ctg	gtg	gaa	tgg	ggc
Q-5	atg	cag	gtg	cag	ctg	gtg	gaa	tgg	ggc
Q-6	atg	cag	gtg	cag	ctg	gtg	gaa	tgg	ggc
E-1	atg	gaa	gtg	cag	ctg	gtg	gaa	tgg	ggc
E-2	atg	gaa	gtg	cag	ctg	gtg	gaa	tgg	ggc
E-3	atg	gaa	gtg	cag	ctg	gtg	gaa	tgg	ggc
E-4	atg	gaa	gtg	cag	ctg	gtg	gaa	agc	ggc
E-5	atg	gaa	gtg	cag	ctg	gtg	gaa	agc	ggc

Fab protein can be synthesized in active form from separate template DNA.

PUREfrex® 2.0 +/- DS supplement / DnaK Mix

↓ + HC and LC (Light Chain) DNA (PCR product) from Herceptin with each ratio
↓ incubation at 37°C or 30°C for 4 h
↓ a), c) SDS-PAGE (0.5 µL of reaction/lane) and Oriole staining
b), d) ELISA

a) SDS-PAGE


Active scFv can be easily synthesized.


PUREfrex® 2.0 +/- DS supplement / DnaK Mix

↓ + scFv DNA (PCR product) from Herceptin
↓ incubation at 37°C for 4 h
↓ a) ↓ spin at 20,000xg for 30 min)
↓ SDS-PAGE (0.25 µL of reaction/lane) and Oriole staining
b) ELISA


c) SDS-PAGE


a) SDS-PAGE


b) SDS-PAGE


b) ELISA

d) ELISA

b) ELISA

AT rich > GC rich @ 5' region
→ High yield of HC

Optimized H/L DNA ratio
→ High yield of active Fab

All tested conditions
→ High yield of active scFv