


Simple and Effective Selection of cyclic peptides by Ribosome display with PURE^{frex}[®]

Kanehisa Kojoh, Kumiko Tsuihiji, Shizue Katoh, Takashi Umeda

GeneFrontier Corporation, Chiba, Japan




Abstract

Cyclic peptides are the peptides internally linked with at least one covalent bond, such as disulfide bond that generate the ring structure. There are several methods reported for the synthesis of cyclic peptide library, but those methods need many steps and complicated techniques (introducing of non-natural amino acids, using of enzymes or chemicals to cyclize peptides etc.), and in most cases, it needs extra work for library screening.

Ribosome display (RD) has many advantages over existing in vitro selection methods, so we developed RD utilizing PURE^{frex}, a fully reconstituted cell-free expression system, and named it as PURE^{frex}RD. Here we report very simple and effective approach for the selection of cyclic peptide using PURE^{frex}RD. Our approach is as follows: A constrained 10mer random peptide library was constructed and screened against a target antigen. The cyclic peptide-ribosome-mRNA complex was formed using the customized PURE^{frex}. After four round selections, following cloning and sequencing resulted in a large variety of the binders with some enrichment. Then, we established screening process by the expression of those peptides in periplasm of *E. coli* as a fusion protein. These fusion proteins had the specific binding affinities for antigen only under the oxidized condition without DTT. This result indicates that the fusion proteins hold the ring structure with disulfide bond, and that the formation of a disulfide bond is necessary for its binding affinity.


As such, we showed that PURE^{frex}RD works for the screening of functional cyclic peptides simply and effectively, and it could be applied for the development of peptide based drugs.

①Best reconstituted cell-free protein synthesis system for *in vitro* selection having the lowest level of RNase contamination: PURE^{frex}



PURE^{frex} was developed as “highly purified PURE system”, and it is clearly more suitable for Ribosome Display. PURE system is a well-known reconstituted in vitro transcription and translation system which consists of purified 36 proteins and *E. coli* ribosomes necessary for transcription, translation and energy recycling. It also contains amino acids, NTPs and *E. coli* tRNA, so the target protein can be synthesized just by addition of the template DNA to the reaction mixture.


②Simple selection with PURE^{frex} RD for cyclic peptide



When the arrest sequence of SecM at 3' terminus is translated in RURE^{frex}, a ribosome can be fixed strongly on mRNA. Also, release factors (RF1, RF2, RF3, RRF) were removed from PURE^{frex}, and oxidized glutathione (GSSG) and disulfide isomerase from *E. coli* (DsbC) in optimized concentration to form disulfide bond were added into PURE^{frex}. As a result, RD complex become highly stable, and cyclized peptide is displayed on ribosome with high efficiency.


③In vitro selection of binders for the target, Erk2

<Condition>
• Target antigen : Biotinylated Erk2 protein




The peptide library with variety of 3×10^{11} was used. Recovered mRNA increased along with the progress of selection round. Final recovery rate against input mRNA after 4th round was about 2.0 %.

④Sequencing after 4th round selection



Some of enrichment were observed and many unique clones were identified from the sequencing (94 clones). In addition, alignment analysis showed that some of the clones were categorized into two groups as in the above. As an overview, arginine and hydrophobic amino acids such as leucine or tryptophan appeared with high frequency.


⑤Simple expression of cyclic peptide as MBP-fusion protein



Some of the enriched clones in sequencing were expressed in periplasm of *E. coli* in very simple construct (FLAG - Cys - 10aa - Cys - c-Myc - His₈), but they could not be detected in western blotting and ELISA (not shown). Therefore, we tried to express it as cyclic peptide-MBP fusion protein as in the above. The data of western blotting (left above) indicated that the expressed fusion proteins were polymerized under oxidized condition in periplasm of *E. coli*. That polymerizations seemed to occur through the SS-bond because polymerization disappeared under reducing condition. In addition, these fusion proteins had the binding activities for the target under the oxidized condition only without DTT (as shown in ELISA data, right above). This result indicates that the fusion proteins hold ring structure with disulfide bond, and that the formation of a disulfide bond is necessary for its binding activity. On the other hand, binder B and D showed very weak signal, which suggest that those clones might be dropped in the screening with this expression & ELISA.


Expression by PURE^{frex}

<Expression reagent>
• PURE^{frex} (customized PURE^{frex})
• DTT : final 2 mM
• Oxidized Glutathione (GSSG) : final 3 mM
• Disulfide isomerase (DsbC) : final 16 μM



Cyclic peptide-MBP fusion proteins were expressed in PURE^{frex}. In this case, the polymerizations were remarkably suppressed than that in *E. coli* expression (as shown in western blotting, left above), and all the binders showed good binding activity (as shown in ELISA, right above). Those results suggest that PURE^{frex} is suitable for the expression of those binders.

⑥Comparison of 1st screening ELISA between PURE^{frex} and *E. coli* expression



Due to the good expression property with PURE^{frex}, it showed remarkably better hit rate than that with *E. coli*. The data suggests that PURE^{frex} is more suitable for the expression and the screening at the initial process, not to miss the potential binders.


⑦Affinity measurement of selected binders by BIACORE

<Condition>
• Sample : Cyclic peptide- MBP-fusion protien (Binder A and C)
• Antigen : Biotinylated Erk2 protein (immobilized RU-162)
• Sensor chip : SA


	k_a (1/Ms)	k_d (1/s)	KD (M)
Binder A	1.48X10 ⁴	5.41X10 ⁻⁴	3.66X10 ⁻⁸
Binder C	1.51X10 ⁴	4.43X10 ⁻⁴	2.93X10 ⁻⁸

Binders were purified by Ni resin. Binding kinetics of binders to Erk2 were determined using surface plasmon resonance with BIACORE. Biotinylated Erk2 was immobilized on a SA sensor chip according to the standard method. At a flow rate of 20 μL/min, five concentrations from 0.3 μM to 4.2 μM were used to record sensorgram, and k_a , k_d and KD were determined by evaluation software.

⑧Application I; For bispecifics (Fab + Cyclic peptide)




⑨Application II; For GFP-fusion protein (Add in to GFP)



⑩Application III; For epitope mapping/mimotope development, etc.

Target antigen : Herceptin (anti-HER2 antibody)

<Supposed application of binder>
• Epitope mapping
• Anti-idiotypic binder (like anti idiotypic antibody)
• Mimotope vaccination



Cyclic peptides were selected against anti-Her2 monoclonal antibody. A lot of binders were obtained from 1st screening ELISA, and some of the enrichments were observed as shown above. The binders didn't bind to control antibody A or B classified into IgG1 or IgG4, respectively. This result shows that those binders bind only to the CDR regions of Herceptin specifically, not to the constant or Fc region.

Summary

1, Simple and effective selection system for the development of cyclic peptide were reported.

• PURE^{frex} (fully reconstituted, customized, the lowest RNase contamination) is a key to success.

2, Actual example of the selection of cyclic peptide against Erk2 as target was performed.

• Many, various, unique binders were obtained in just 4 round panning.
• 10⁻⁸ KD binders were selected without affinity maturation.

3, Simple screening is possible with PURE^{frex}.

• Only with natural amino acids, expression is easy.
• Many clones can be expressed for the screening in parallel.
• PURE^{frex} makes it more effective comparing to *E. coli*.

4, Unique applications were proposed.

• Cyclic peptide can add another binding activity to the existing biologics to make it bispecific in very simple manner.
• Targeting existing antibodies, epitope mapping, mimotope development, anti-idiotypic binder (like anti-idiotypic antibody) development can be done with PURE^{frex} RD.

(On-going)

• Unique cyclic peptide for PD-1 (and for other therapeutic targets) are under validation.

• DRP (Disulfide Rich Peptide) screening is on-going aiming for the development of orally available biologics.

<Acknowledgments>
We thank Mr. Keita Iguchi in KANEKA corporation for BIACORE analysis.

Contact
URL: www.genefrontier.com
E-mail: contactus@genefrontier.com