


Selection and affinity maturation of cyclic peptide against CTLA-4 with PURE_{frex}RD

Kumiko Tsuihiji; Shizue Katoh; Mikiko Nakamura; Satoshi Murakami; Kanehisa Kojoh

GeneFrontier

GeneFrontier Corporation, Chiba, Japan (<https://www.genefrontier.com/en/>)

15th Annual PEGS BOSTON (April 8-12, 2019)



Abstract

[Background]

Ribosome display (RD) is the simplest method among existing *in vitro* selection methods. We developed RD by optimizing the composition of PURE_{frex}, for instance, removing release factors from the system, adjusting the concentration of metal ions and ribosome, etc. We named it as PURE_{frex}RD. We have already reported that a lot of cyclic peptides and antibodies were successfully selected simply and effectively by PURE_{frex}RD. This time, we carried out *in vitro* selection using PURE_{frex}RD to select cyclic peptides binding to cytotoxic T-lymphocyte-associated protein 4 (CTLA-4) with inhibitory activity to the interaction between CTLA-4 and CD80, which can be one of the new modality as an immune check point inhibitor.


[Methods and Results]

Cyclic peptide library having 12mer random sequences (Cys - X12 - Cys) was screened against a biotinylated CTLA-4-Fc fusion protein. After 3 rounds of selection, enriched mRNAs were subcloned, and a lot of enriched clones with variety were obtained from sequence analysis. The enriched clones were examined for their binding activity to CTLA-4 by RD pull down assay, and 3-02 clone showed specific binding activity to CTLA-4. To evaluate the affinity to CTLA-4 and the inhibitory activity to the interaction between CTLA-4 and CD80, 3-02 clone was expressed as MBP-fusion protein in *E. coli* and purified. Purified 3-02 clone showed very low affinity and very low inhibitory activity, so affinity maturation of 3-02 clone by off-rate selection using PURE_{frex}RD was performed to improve its affinity and inhibitory activity. Mutated library against 3-02 clone was prepared by error prone PCR. After 3 rounds of off-rate selection over several hours (1st round; 2 hours, 2nd round; 19 hours, 3rd round; 67 hours), mutated clones were expressed in *E. coli*, and their binding activities were evaluated by ELISA. The purified clones showed the higher activity in ELISA, and each EC₅₀ was determined for comparison in specific activity of each clone to the original 3-02 clone. Finally, we have obtained several high affinity clones by approximately 50-fold more than the original 3-02 clone.

[Conclusion]

PURE_{frex}RD will be applied to simple and effective *in vitro* selection of cyclic peptide and will contribute to the development of cyclic peptide based drugs.

In vitro selection by PURE_{frex}RD




Cyclic peptide library having 12mer random sequences (Cys - X12 - Cys) was screened against biotinylated CTLA-4-Fc protein. Recovered mRNA increased along with the progress of selection round.

Sequencing after 3rd round selection

Number of conc.	clone name	Sequence														
		Cys	1	2	3	4	5	6	7	8	9	10	11	12	Cys	
51	3-02	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
6	3-51	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
4	3-33	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
3	3-80	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
2	2-23	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
2	3-90	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
2	3-38	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y

Some of enrichment were observed from the sequencing (94 clones). Each color represents specific amino acid.

Confirmation of binding activity by RD pull down assay



Only 3-02 binder which was the most enriched clone bound to the CTLA-4 region in CTLA-4-Fc protein specifically (red-square), and another clones bound to the Fc region.

Affinity measurement of 3-02 binder by BLItz system

<Instrument for measurement>
• BLItz system (FORTE Bio)


<Condition>

- Sample : Cyclic peptide- MBP-fusion protein
- Antigen : Biotinylated CTLA-4-Fc protein
- Biosensor : SA

	k_a (1/Ms)	k_d (1/s)	KD (M)
3-02	3.18×10^3	1.23×10^{-2}	3.88×10^{-6}


3-02 binder was expressed as MBP-fusion protein in *E. coli* and purified by affinity resin. Binding kinetics to CTLA-4-Fc was determined using BLItz system. Biotinylated CTLA-4-Fc was immobilized on a SA Biosensor according to the standard method.

Affinity maturation by PURE_{frex}RD (Off-rate selection)



In off-rate selection, RD complexes bound to CTLA-4-immobilized streptavidin beads were washed with large excess of free CTLA-4 to prevent the rebinding of RD complexes to the beads. The higher-affinity binders can be retained on the beads in a longer washing time. Library was prepared by error prone PCR using 3-02 clone as template DNA.

Evaluation of amino acid sequence and ELISA in MBP-fusion format after affinity maturation by PURE_{frex}RD



It seems that mutating of position 6 residue was very important to have the higher binding activity, and position 2,3,5,9,11,12 residues would be essential for specific binding to CTLA-4.


Measurement of affinity (EC₅₀) as MBP-fusion peptide

E. coli expression host : BL21(DE3)

Culture condition : 30


Induction : final 0.1 mM IPTG

Purification : His tag purification





The high affinity clones (27 clones) from ELISA were purified, and EC₅₀ of the each were measured. As a result, all clones except No.27 showed the higher EC₅₀ than original. In particular, No.19, No. 24 (red square) showed about 50 times higher affinities than the original.

Cell based binding assay with the synthetic peptide (FACS)



No.19 binder as synthetic peptide (N-term biotinylated) bound to CTLA-4 on the cells.

Measurement of inhibition (IC₅₀)




No.19 binder showed the inhibitory activity against interaction between CD80 and CTLA-4 in both MBP-fusion and fully synthetic peptide format, but original clone (3-02) showed very low activity. However, the inhibitory activity of the synthetic peptide was lower than MBP-fusion. It seemed that the activity of synthetic peptide was influenced by the characters of peptide such as solubility, isoelectric point and conformational fluctuation.

Summary

- 1, Functional cyclic peptides against CTLA-4 were selected easily and rapidly by PURE_{frex}RD.
- 2, In affinity maturation of CTLA-4 binders, we succeeded to select lead cyclic peptide with good EC₅₀/IC₅₀.

E-mail: contactus01@genefrontier.com

Best reconstituted cell-free protein synthesis system for *in vitro* selection having the lowest level of RNase contamination: PURE_{frex}




- Points
- Fully reconstituted.
 - No cells, No bias.
 - Fast, Easy, Pure.
 - High throughput.
 - Lower contaminants.
- Very low level of RNase, patented as JP4931135

Best suited for Ribosome display and mRNA display

PURE_{frex} is the reconstituted *in vitro* transcription and translation system which consists of purified 36 proteins and *E. coli* ribosomes necessary for transcription, translation and energy recycling. It also contains amino acids, NTPs and *E. coli* tRNA, so the target protein can be synthesized just by addition of the template DNA to the reaction mixture.


Simple selection with PURE_{frex}RD for cyclic peptide

Ribosome Display complex with cyclic peptide



When the arrest sequence of SecM at 3' terminus is translated in PURE_{frex}, a ribosome can be fixed strongly on mRNA. Also, release factors (RF1, RF2, RF3, RRF) were removed from PURE_{frex}, and oxidized glutathione (GSSG) and disulfide isomerase from *E. coli* (DsbC) in optimized concentration to form disulfide bond were added into PURE_{frex}. As a result, RD complex become highly stable, and cyclized peptide is displayed on ribosome with high efficiency.

Immune checkpoint via CTLA-4



When CTLA-4 bind to B7, T-cells are led to anergy. Anergic T-cells have limited effector function. Anti-CTLA-4 neutralizing cyclic peptides inhibit CTLA-4 binding to B7 and promote T-cell activation.