

In vitro selection and affinity maturation of CTLA-4 binding cyclic peptide with PUREfrexRD, and conversion of the peptide to small molecules.

PUREfrexRDによるCTLA-4結合環状ペプチドの試験管内選択と親和性向上および環状ペプチドの低分子化合物への変換



GeneFrontier

對比地 久美子¹, 加藤 静恵¹, 中村 美紀子¹, 村上 智史¹, 高島 一², 河西政次², 毛利 睦², 吉森篤史³, 〇古城 周久¹
¹ジーンフロンティア (株) ² (株) PRISM BioLab, ³理論創薬研究所

Abstract

[Background]


In the drug development, small molecule drugs have a lot of advantage over peptides in terms of oral administration and intracellular migration. However, since diversity of small molecule library and selection methods are limited, it is difficult to get the promising drugs that has the inhibitory activity for protein-protein interactions. On the other hands, we have reported that the functional cyclic peptides can be easily selected / screened by the ribosome display method (PUREfrexRD) using PUREfrex. Here, we report selection and affinity maturation of functional cyclic peptide that inhibited the binding between CD80 and CTLA-4, and conversion of the cyclic peptide to small molecules.

[Method and Results]

A lot of cyclic peptides that specifically binds to CTLA-4 was obtained using PUREfrexRD (3 rounds selection) from cyclic peptide library consisting of a random sequence of 12 residues. Next, a cyclic peptide having inhibitory activity was selected and applied for affinity maturation (off-rate selection) by PUREfrexRD. As a result, a high affinity cyclic peptide having an improved affinity of about 630 times was obtained. Furthermore, the comprehensive information related to the correlation between amino acid sequence and binding activity was obtained from a series of experiments such as identification of essential amino acids for binding by alanine scanning, exhaustive expression of recombinant mutants and ELISA. Based on these data, the 3D-structure of the cyclic peptide was predicted in silico, and about 350 small molecules similar to the simulated cyclic peptide structure were picked out from PRISM BioLab's small molecule library, and binding inhibition activity were estimated by ELISA. As a result, a several small molecules having inhibitory activity could be obtained.

It indicates that PUREfrexRD is an effective way for obtaining functional cyclic peptides and their structure-activity relationship information, and conversion from cyclic peptides to small molecules.

In vitro selection by PUREfrexRD



Cyclic peptide library having 12mer random sequences (Cys - X12 - Cys) was screened against biotinylated CTLA-4-Fc protein. Recovered mRNA increased along with the progress of selection round.

Sequencing after 3rd round selection

Number of conc.	clone name	Sequence													
		Cys	1	2	3	4	5	6	7	8	9	10	11	12	Cys
51	3-02	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
6	3-51	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
4	3-33	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
3	3-80	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
2	2-23	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
2	3-90	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
2	3-38	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y


Some of enrichment were observed from the sequencing (94 clones). Each color represents specific amino acid.

Measurement of affinity (EC50) as MBP-fusion peptide after affinity maturation (off-rate selection) by PUREfrexRD

(Off-rate selection)
 Library: Error prone mutants
 Wash time: 1st-2H, 2nd-19H, 3rd-67H

(Expression)
E. coli host : BL21(DE3)
 Culture condition : 30

(Purification)
 His tag purification




Original	1	2	3	4	5	6	7	8	9	10	11	12	EC50 µM
No.19													0.075
No.10													0.002
No.24													0.005
No.20													0.180
No.3													0.180
No.18													0.225
No.16													0.217
No.6													0.222
No.4													0.225
No.11													0.268
No.5													0.300
No.21													0.322
No.14													0.342
No.7													0.352
No.2													0.369
No.8													0.414
No.6													0.508
No.11													0.500
No.23													0.674
No.14													0.703
No.17													0.808
No.13													0.832
No.9													0.867
No.12													0.869
No.1													1.005
No.22													1.403
No.16													1.603
No.27													4.181

The high affinity clones (27 clones) from ELISA were purified, and EC50 of the each were measured. As a result, all clones except No.27 showed higher EC50 than the original (3-02). In particular, No.19, 10, 24 (red square) showed about 50 times higher affinities than the original.


Identification of amino acids required for binding by alanine scanning

(Construct)
 MBP-fusion peptide




Amino acids of position 2, 3, 5, 6, 9, 11, 12 (red circle) were necessary for binding to CTLA-4.

Cell based binding assay with the synthetic peptide (FACS)




No.19 binder as synthetic peptide (N-term biotinylated) bound to CTLA-4 on the cells.

Measurement of inhibition (IC50) by inhibitory ELISA





No.19 binder showed the inhibitory activity against interaction between CD80 and CTLA-4 in both MBP-fusion and fully synthetic peptide format. The inhibitory activity of the synthetic peptide was lower than MBP-fusion.

Designing small molecules based on cyclic peptides




J.Med.Chem., 45, 1396 (2002)



Small molecules mimicking beta-turn derived from cyclic peptides


Three amino acids are introduced into the R position (red square) of the beta-turn peptide-mimetic scaffold in order from the N-terminus of the cyclic peptide.

Inhibitor Screening assay and Measurement of IC50



Several small molecules with inhibition activity were obtained. IC50 of Pm2 with the highest activity showed about 1/20 of the synthetic peptide.

Comparison of affinity among Pm2 analogs




Pm2 was designed from 3 amino acids including essential 2 amino acids for binding at C-terminal of peptide. Two analogs with adduct or substitution of 1 amino acid showed the decrease in IC50.

Summary

- 1, Functional cyclic peptides against CTLA-4 were selected easily and rapidly by PUREfrexRD.
- 2, In affinity maturation of CTLA-4 binder, we succeeded to select lead cyclic peptide with good EC50/IC50.
- 3, Essential amino acids for binding were identified by alanine scanning.
- 4, Based on structure-activity relationship information, peptide-mimetic small molecules were designed, and successfully selected with good inhibitory activity.

E-mail: contactus01@genefrontier.com

Best reconstituted cell-free protein synthesis system for in vitro selection having the lowest level of RNase contamination: PUREfrex



- Points
- Fully reconstituted.
 - No cells, No bias.
 - Fast, Easy, Pure.
 - High throughput.
 - Lower contaminants.
- Very low level of RNase, patented as JP4931135


Best suited for Ribosome display and mRNA display

Shimizu Y et al., Nat. Biotechnol. (2001) vol 19, p751-755

PUREfrex is the reconstituted in vitro transcription and translation system which consists of purified 36 proteins and E. coli ribosomes necessary for transcription, translation and energy recycling. It also contains amino acids, NTPs and E. coli tRNA, so the target protein can be synthesized just by addition of the template DNA to the reaction mixture.


Simple selection with PUREfrexRD for cyclic peptide

Ribosome Display complex with cyclic peptide



When the arrest sequence of SecM at 3' terminus is translated in PUREfrex, a ribosome can be fixed strongly on mRNA. Also, release factors (RF1, RF2, RF3, RRF) are removed from PUREfrex, and oxidized glutathione (GSSG) and disulfide isomerase from E. coli (DsbC) in optimized concentration to form disulfide bond are added into PUREfrex. As a result, RD complex become highly stable, and cyclized peptide is displayed on ribosome with high efficiency.

Immune checkpoint via CTLA-4



When CTLA-4 bind to B7, T-cells are led to anergy. Anergic T-cells have limited effector function. Anti-CTLA-4 neutralizing cyclic peptides inhibit CTLA-4 binding to B7 and promote T-cell activation.